

LABSTAC

STRONG ACID AND ALKALI STORAGE CABINET CBT21-45L



labstac.com

STRONG ACID AND ALKALI STORAGE CABINET

CBT21-45L

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT21-45L STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

Model	CBT21-45L
Capacity	45 gallon (170L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1650mm
Package Size	1175x545x1785 mm
Gross Weight	72 Kg



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com