

STRONG ACID AND ALKALI STORAGE CABINET CBT21-4L



STRONG ACID AND ALKALI STORAGE CABINET

CBT21-4L

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT21-4L STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

Model	CBT21-4L
Capacity	4 gallon (15L)
Material	pp
Adjustable Shelf	Optional
Door	Single / Manual
External Size	430x430x500 mm
Package Size	515x515x685 mm
Gross Weight	16 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com