

**LABSTAC**

# **STRONG ACID AND ALKALI STORAGE CABINET CBT21-90L**



labstac.com

# STRONG ACID AND ALKALI STORAGE CABINET

## CBT21-90L

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## CBT21-90L STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



### SPECIFICATIONS

Model	CBT21-90L
Capacity	90 gallon (340L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x860x1650 mm
Package Size	1173x943x1785 mm
Gross Weight	107 Kg

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)