

# HEMATOCRIT CENTRIFUGE CEN34-1200HE



# HEMATOCRIT CENTRIFUGE CEN34-1200HE

These laboratory multi purpose Medical Centrifuges are ideal solution for high-end research experiment such as PRP/PPP, cell washing, HCT, blood type test etc. in the biology, radio-immunity, medical and other research fields. These versatile centrifuges have precise temperature control with highly efficient double loop control design which allows fast heating and cooling.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

## CEN34-1200HE HEMATOCRIT CENTRIFUGE

LED display

Microprocessor controller

Brushless DC motor

Electronic door interlock

### SPECIFICATIONS

Model	CEN34-1200HE
Maximum Capacity (No of tubes x Vol.)	24 Capillary Tubes
Maximum Speed	12000 rpm
Speed Accuracy	± 20 rpm
Maximum RCF	15800xg
Time Range	1~99 min 59 s : 59 s
Drive System	Direct, brushless motor drive
Overall Dimension	10.6x14.2x8.3 / 27x36x21 mm
Noise Level	≤ 55 dB(A)
Net. Weight	14 kg
Max. Power Consumption	500 W
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)