

# BLOOD TYPING CENTRIFUGE CEN35-0384BT



# BLOOD TYPING CENTRIFUGE CEN35-0384BT

These laboratory multi purpose Medical Centrifuges are ideal solution for high-end research experiment such as PRP/PPP, cell washing, HCT, blood type test etc. in the biology, radio-immunity, medical and other research fields. These versatile centrifuges have precise temperature control with highly efficient double loop control design which allows fast heating and cooling.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

## CEN35-0384BT BLOOD TYPING CENTRIFUGE

- Brushless DC motor
- LED display
- Microprocessor controller
- The inner chamber is made up of stainless steel
- Electronic door interlock
- Faults are diagnosed automatically
- Operating parameters can be modified during the operation
- Stable operation



## SPECIFICATIONS

Model	CEN35-0384BT
Maximum Capacity (No of tubes x Vol.)	24 gel card
Maximum Speed	3840 rpm
Maximum RCF	1790xg
Drive System	Direct, induction motor drive
Overall Dimension	14.8x14.2x11.8 in / 37.5x36x30 mm
Noise Level	≤ 60 dB(A)
Weight (Net/Gross)	21 kg
Max. Power Consumption	100 W
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Rotor Capacity
2201310005	No. 1 Swing-out Rotor	12 gel cards
2201310006	No. 2 Swing-out Rotor	24 gel cards

***LABSTAC***

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)