

# CELL-WASHING CENTRIFUGE CEN35-0470CW



# CELL-WASHING CENTRIFUGE CEN35-0470CW

These laboratory multi purpose Medical Centrifuges are ideal solution for high-end research experiment such as PRP/PPP, cell washing, HCT, blood type test etc. in the biology, radio-immunity, medical and other research fields. These versatile centrifuges have precise temperature control with highly efficient double loop control design which allows fast heating and cooling.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

## CEN35-0470CW CELL-WASHING CENTRIFUGE

LED display

Microprocessor controller

Brushless DC motor

Electronic door interlock

### SPECIFICATIONS

Model	CEN35-0470CW
Maximum Capacity (No of tubes x Vol.)	12x7 ml (SERO Rotor)
Maximum Speed	4700 rpm
Speed Accuracy	± 20 rpm
Maximum RCF	2000xg
Drive System	Direct, induction motor drive
Overall Dimension	14.8x14.2x11.8 in / 37.5x36x30 mm
Noise Level	≤ 55 dB(A)
Weight (Net/Gross)	16 kg
Max. Power Consumption	200 W
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz

### OPTIONAL ACCESSORIES

Accessory Code	Name	Rotor Capacity
2201311005	No. 1 HLA Rotor/WBC Rotor	12x1.5/2.0 ml
2201311006	No. 2 SERO Rotor/RBC Rotor	12x7 ml

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)