

# MICROCENTRIFUGE CEN43-060



# MICROCENTRIFUGE CEN43-060

This compact, powerful micro centrifuge is useful to separate or sediment small volume samples at high speed. The short-spin key, soft brake, FastTemp features optimizes the overall function of the centrifuges. Faults are diagnosed automatically and displays fault code. These hybrid micro-centrifuge can run PCR tubes, Capillary tubes along with standard tubes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

## CEN43-060 MICROCENTRIFUGE

Quick release of rotor for easy rotor exchange

Gentle braking system for quick sample processing



## SPECIFICATIONS

Model	CEN43-060
Maximum Capacity (No of tubes x Vol.)	8x1.5 ml / 2.0 ml
Speed Range	4000 & 6000
Maximum RCF	1200xg & 2000xg
Overall Dimension	135x147x120 mm
Weight (Net/Gross)	0.7 kg
Power Supply	AC100-240V

## OPTIONAL ACCESSORIES

Accessory Code	Name	Capacity
2202614005	Fixed-Angle Rotor	8x1.5/2.0ml
2202614006	Strip Rotor	4-8x0.2ml
2202614007	Adapter	0.5-0.6ml tube
2202614008	Adapter	0.2ml tube



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)