

DIFFERENTIAL SCANNING CALORIMETER DSC13-600B



DIFFERENTIAL SCANNING CALORIMETER

DSC13-600B

DSC13-600B DIFFERENTIAL SCANNING CALORIMETER

Industrial level 7-inch touch screen, display information rich

New furnace body structure, can be matched with a variety of refrigeration devices

USB communication interface, strong versatility, reliable communication, support self recovery connection function

Automatic switching of two-way atmosphere flow, fast switching speed and short stable time. At the same time, a protective gas input is added

The software is simple and easy to operate

SPECIFICATIONS

Model	DSC13-600B
Temperature range	-150 °C ~ 600 °C
Temperature resolution	0.01°C
Temperature fluctuation	± 0.01°C
Temperature repeatability	± 0.1°C
Heating rate	0.1~80°C/min
Constant temperature time	0~400
DSC range	0~±600mW
DSC resolution	0.001 uW
DSC accuracy	0.001 mW
Working power supply	Ac220 v 50 Hz
Atmosphere control gas	nitrogen, oxygen
Cooling equipment	Nitrogen cooling
Gas flow	0-300 mL/min
Gas pressure	0.2 MPa
Display mode	24 bit color, 7 inch LCD touch screen display
Data interface	standard USB interface

The logo for Labstac, featuring the word "LABSTAC" in a bold, sans-serif font. The letters "LAB" are red, and "STAC" is black.

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com