

TOP PRESS FREEZE DRYER FDR21-01LTE



TOP PRESS FREEZE DRYER FDR21-01LTE

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

FDR21-01LTE TOP PRESS FREEZE DRYER



SPECIFICATIONS

Model	FDR21-01LTE
Type	Top Press with Manifold type electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	0.1 m ²
Bulk capacity	1.0 L
Tray	φ180mm, 4pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com