

# TOP PRESS FREEZE DRYER FDR21-01TM



# TOP PRESS FREEZE DRYER FDR21-01TM

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## FDR21-01TM TOP PRESS FREEZE DRYER



### SPECIFICATIONS

Model	FDR21-01TM
Type	Top Press with Manifold Type
Condenser Temperature	-110°C
Freeze Drying area	0.1 m <sup>2</sup>
Bulk capacity	1.0 L
Tray	φ180mm, 4pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)