

LABSTAC

T-MANIFOLD FREEZE DRYER FDR21-18LMT



labstac.com

T-MANIFOLD FREEZE DRYER FDR21-18LMT

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.

Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

FDR21-18LMT T-MANIFOLD FREEZE DRYER



SPECIFICATIONS

Model	FDR21-18LMT
Type	Manifold type Electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	0.18 m ²
Bulk capacity	1.8 L
Tray	φ200mm, 6pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com