

# T-MANIFOLD FREEZE DRYER FDR21-27LMT



# T-MANIFOLD FREEZE DRYER FDR21-27LMT

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.

Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## FDR21-27LMT T-MANIFOLD FREEZE DRYER



### SPECIFICATIONS

Model	FDR21-27LMT
Type	Manifold Type
Condenser Temperature	-110°C
Freeze Drying area	0.18 or 0.27 m <sup>2</sup>
Bulk capacity	1.8 or 2.7 L
Tray	φ240mm, 4 or 6pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)