

TOP PRESS FREEZE DRYER FDR21-27LST



TOP PRESS FREEZE DRYER FDR21-27LST

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

FDR21-27LST TOP PRESS FREEZE DRYER



SPECIFICATIONS

| | |
|-----------------------|-----------------------------|
| Model | FDR21-27LST |
| Type | Standard type |
| Condenser Temperature | -110°C |
| Freeze Drying area | 0.18 or 0.27 m ² |
| Bulk capacity | 1.8 or 2.7 L |
| Tray | φ240mm, 4 or 6pcs |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com