

# STANDARD FREEZE DRYER FDR21-36STE



# STANDARD FREEZE DRYER FDR21-36STE

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## FDR21-36STE STANDARD FREEZE DRYER



### SPECIFICATIONS

Model	FDR21-36STE
Type	Standard type electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	: 0.36 m <sup>2</sup>
Bulk capacity	3.6 L
Tray	φ240mm, 8pcs

# **LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)