

IN SITU FREEZE DRYER FDR22-03M



IN SITU FREEZE DRYER FDR22-03M

Compact, stainless steel unit with casters for easy movement, food type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.

FDR22-03M IN SITU FREEZE DRYER



SPECIFICATIONS

Model	FDR22-03M
Type	Food freeze Dryer
Lowest condenser temperature	-50 °C
Drying area	0.3m ²
Tray temperature	heating to +50 °C
Trays load capacity	3L

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com