

# IN SITU FREEZE DRYER FDR22-56M



# IN SITU FREEZE DRYER FDR22-56M

Compact, stainless steel unit with casters for easy movement, food type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.

## FDR22-56M IN SITU FREEZE DRYER



### SPECIFICATIONS

Model	FDR22-56M
Type	Food freeze Dryer
Lowest condenser temperature	-50 °C
Drying area	0.56m <sup>2</sup>
Tray temperature	heating to +50 °C
Trays load capacity	6L

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)