

# IN SITU FREEZE DRYER FDR23-02M



# IN SITU FREEZE DRYER FDR23-02M

Silicone Oil Heating In-situ pilot freeze Dryer changes the traditional complicated operation, prevents the pollution of the material realizes the dry sublimation automation. It is suitable for unit production

Used in Food processing, Pharmaceutical, and Medical sector.

Also known as Silicone Oil Heating Lyophilizer.

## FDR23-02M IN SITU FREEZE DRYER

One-piece design, compact equipment installation

The maintenance orientation is concentrated on one side, which is convenient for maintenance

A variety of security interlocking functions to fully protect the product

The equipment is conveniently connected to the clean room, and can be connected to the isolator

Realize dual control of clean room and maintenance room, but interlocking with each other

Simple operation, all run automatically



## SPECIFICATIONS

Model	FDR23-02M
Type	Pilot In-situ freeze Dryer
Condenser Temperature	-75°C
Drying area	2m <sup>2</sup>
Trays load capacity	30kg
Liquid loading capacity	30L
Computer connection	RS232

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)