

**LABSTAC**

# STANDARD FUME HOOD FUH12-1500



labstac.com

# STANDARD FUME HOOD FUH12-1500

Created to provide protection from toxic fumes and hazardous volatile vapours, and are reliable and easy to use. Inbuilt ventilation system lowers the risk of fire and explosion. For metal analysis, pharmacy and research field along with safety measures, this product becomes unique choice for the buyers to purchase according to their needs.

Used in Metal analysis, Pharmacy, Research Field, Commercial, Medical, Laboratory.

Also known as Standard Fume Hood, Laboratory Fume Hood, Chemical Fume Hood.

## FUH12-1500 STANDARD FUME HOOD

Interior body includes melamine board with acid and alkali resistance function

LCD display

Exterior body is of cold-roll steel with anti-bacterial powder

Working chamber is made of solid chemical-resistant material

UV lamp to sterilize the working area

Adjustable air velocity with electrical controlled glass

Adjustable blower speed of a centrifugal fan



## SPECIFICATIONS

Model	FUH12-1500
Air Velocity	0.3~0.8 m/s
Max Opening	520 mm
Work Surface Height	750 mm
System Exhaust Volume	990m <sup>3</sup> /h
Front Window	Motorized, 5 mm toughened glass; Anti-UV Exterior Body: cold rolled steel with anti-bacteria powder coating.
Work Table Material	Chemical resistant Phenolic resin
LED Lamp	16Wx 1
UV Lamp	20 W Emission of 253.7 nanometers for most efficient decontamination
Blower	Built-in centrifugal blower, Speed adjustable with 9 levels
Standard Accessory	Water tap, Gas tap, Water sink, base stand, UV lampx2, LED lamp, Socket, 4.3 meters exhaust duct, Active carbon filter, foot switch, Waterproof Socket- 2 pcs
Optional Accessory	HEPA filter
Internal Size (WxDxH)	1380x730x745 mm
External Size (WxDxH)	1500x840x2150 mm
Package (Main Body)	1650x1120x1680 mm
Noise Level	≤60 dB
Gross Weight	290 kg
Consumption	500 W
Power Supply	AC220V±10%, 50/60Hz 110V+10%, 60Hz

## ACCESSORIES

Accessory Code	Name
5900509005	Fluorescent lamp
5900509006	UV lamp
5900509007	Base stand
5900509008	Waterproof socket
5900505005	water tap

Accessory Code	Name
5900505006	Gas tap
5900505007	Water sink
5900509009	4 meters PVC exhaust duct
5900505011	Active carbon filter
5900509010	Foot switch

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)