

AUTOMATIC AIR GENERATOR GEN11-2000



AUTOMATIC AIR GENERATOR GEN11-2000

Air Generator is designed for removing hydrocarbons from the air stream and produces pure, hydrocarbon free and dry air, with a very low dew point. Long service life, easy maintenance, abrasion resistant, cost effective, performance, quick and easy servicing are some of the features which makes it convenient for the laboratory.

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

GEN11-2000 AUTOMATIC AIR GENERATOR

High precision output pressure

Water automatically drain out

Two levels stainless steel made purifications

Large inside stainless steel air tank

Have two levels pressure stable valve. Output air pressure much stable

Have over pressure protection, low voltage start up, DC, greatly protect inside compressor



SPECIFICATIONS

Model	GEN11-2000
Output Pressure	0-0.4MPa
Noise	<35dB
Dimension	420x260x350 mm
Net Weight	22 kg
Consumption Power	150 W
Output Flow Rate	0-2000 ml/min
Environmental Temperature	0-45°C
Humidity	≤85%

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com