

AUTOMATIC NITROGEN GENERATOR GEN32-1000



AUTOMATIC NITROGEN GENERATOR

GEN32-1000

Designed in such a way that not only removes the hassle of having to change over cylinders, it provides an uninterrupted supply of gas at a consistent purity. Features like low operating cost, environment friendly, safe, energy efficient makes it an ideal choice for laboratory.

Used in Analytical chemistry, Spectroscopy, Pharmaceutical, General laboratory, Food Packaging, Reducing fire hazards, Manufacturing of stainless steel, Wine bottling, LC-MS Gas Chromatography, Fire protection, Glass Industry..

Also known as Laboratory Nitrogen Generator.

GEN32-1000 AUTOMATIC NITROGEN GENERATOR

Automatic control, constant pressure, constant flow

Ultra high molecular membrane separation technique, PCAN carrier and precious metal catalyst

Deoxidation and Dehydration device, high purity

Compact, frees up valuable space

Easy operation, eliminate the need for expensive gas cylinder

Safe and reliable



SPECIFICATIONS

Model	GEN32-1000
Nitrogen purity	Oxygen content < 3PPM
Flow Rate	0-1000 ml / min
Output Pressure	0-0.4 Mpa
Pressure Stability	< 0.003 Mpa
Dimension	460x360x350 mm
Net Weight	20 kg
Consumption Power	150 W
Power Supply	220 V ± 10%, 50 Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com