

# AUTOMATIC HYDROGEN AIR GENERATOR GEN41-0500



# AUTOMATIC HYDROGEN AIR GENERATOR

## GEN41-0500

Hydrogen Air Generator provides hydrogen and pure air simultaneously. They utilize PEM technology to produce high-purity hydrogen gas, no additives required for corrosion and pollution free operation, making our products an ideal choice for laboratory.

Used in Hydrogen & air supply for fuel cells, Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

### GEN41-0500 AUTOMATIC HYDROGEN AIR GENERATOR

- Stainless steel air tank inside
- Deoxidization device
- Output gas purity is higher
- Inside machine use silicon rubber ring ( lower content sulfur)
- Increased gas output quality
- Compact design, easy using, automatically operation, safety



### SPECIFICATIONS

Model	GEN41-0500
Hydrogen Flow Rate	0-500 ml/min
Air Flow Rate	0-2000 ml/min
Hydrogen Pressure	0.4 Mpa
Dimension	470x260x380 mm
Net Weight	28 kg
Consumption Power	400 W



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)