

AUTOMATIC HYDROGEN AIR GENERATOR GEN41-0500



AUTOMATIC HYDROGEN AIR GENERATOR

GEN41-0500

Hydrogen Air Generator provides hydrogen and pure air simultaneously. They utilize PEM technology to produce high-purity hydrogen gas, no additives required for corrosion and pollution free operation, making our products an ideal choice for laboratory.

Used in Hydrogen & air supply for fuel cells, Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

GEN41-0500 AUTOMATIC HYDROGEN AIR GENERATOR

Stainless steel air tank inside

Deoxidization device

Output gas purity is higher

Inside machine use silicon rubber ring (lower content sulfur)

Increased gas output quality

Compact design, easy using, automatically operation, safety



SPECIFICATIONS

Model	GEN41-0500
Hydrogen Flow Rate	0-500 ml/min
Air Flow Rate	0-2000 ml/min
Hydrogen Pressure	0.4 Mpa
Dimension	470x260x380 mm
Net Weight	28 kg
Consumption Power	400 W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com