

THREE FUNCTION ELECTRIC BED HBD14-03



THREE FUNCTION ELECTRIC BED HBD14-03

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD14-03 THREE FUNCTION ELECTRIC BED

Steel Framework, top class ABS & PP plastic

3 functions controlled by remote control

Castor with central brake system, wheel diameter: 150mm



SPECIFICATIONS

Model	HBD14-03
Function	3
Back-rest lifting angle	0° - 70°
Knee Gatch	0°-30°
HI – LOW	665mm – 450 mm
Castor	360° swivel castors (Dia. 150 mm), with central brake under the foot board.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1940x900 mm
Safe load capacity	< 250 kg
Standard accessory	IV pole
Dimension	2160Lx1120Wx440-655 mm
Net Weight	139 kg
Patient surface	1935Lx900W mm
Height adjustment	440-655 mm
IV pole	Yes

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@m.labstac.com | Website: labstac.com