

# MULTI- FUNCTION ELECTRIC BED HBD16-M



# MULTI- FUNCTION ELECTRIC BED HBD16-M

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD16-M MULTI- FUNCTION ELECTRIC BED

Steel Framework, top class ABS & PP plastic

Functions activated by electric controlling adopt luxurious membrane switch, easily operated and pressed comfortably

With the "reset" button to make the bed be flat from any position

A lock button with indicator light to show work state

Protective railing with angle indicators

castor with central brake system, wheel diameter: 150mm



## SPECIFICATIONS

Model	HBD16-M
Function	Multi
Back-rest	0° - 65°
Knee Gatch	0° - 40°
Trendelenburg	0 - 12°
Reverse trendelenburg	0 – 12°
HI – LOW	745 – 477 mm
Castor	360° swivel castors (Dia. 150 mm), with central brake.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1940x900 mm
Safe load capacity	< 250 kg
Standard accessory	IV pole
Dimension	2100Lx1050Wx485-755 mm
Net Weight	139 kg
Height adjustment	485-755 mm
IV pole	Yes



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)