

MANUAL BED FOUR CRANK 5 FUNCTION HBD27-05



MANUAL BED FOUR CRANK 5 FUNCTION

HBD27-05

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD27-05 MANUAL BED FOUR CRANK 5 FUNCTION

ABS head & foot boards can be removable.

Aluminium alloy side rail, 6-staff, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish.

Bed framework made from mild steel ERW rectangular tubes.

Covered castors dia125mm with brake

With telescopic SS IV pole and prevision

SPECIFICATIONS

Model	HBD27-05
Function	5
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Height adjustment	500-750 mm (without mattress)
IV pole	Yes
Dimension	2080x970x500 mm

The logo for Labstac, featuring the word "LABSTAC" in a bold, sans-serif font. The letters "LAB" are in red and "STAC" are in black.

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com