

FLAT MEDICAL BED HBD2B-FM



FLAT MEDICAL BED HBD2B-FM

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD2B-FM FLAT MEDICAL BED

ABS head & foot boards can be removable. with bumper, 8kg/pair, color :blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

4 pcs Individual lockable Taiwan Secure castors at 4 corners, diam 5"

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD2B-FM
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
IV pole	Yes
Dimension	2160Lx980Wx500H mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com