

MANUAL 2

FUNCTION MEDICAL

BED HBD2E-02



MANUAL 2 FUNCTION MEDICAL BED HBD2E-02

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD2E-02 MANUAL 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

| | |
|-------------------------|--------------------|
| Model | HBD2E-02 |
| Function | 2 |
| Back-rest lifting angle | 0°-80°±5° |
| Foot-rest lifting angle | 0°-40°±5° |
| IV pole | Yes |
| Dimension | 2150Lx980Wx500H mm |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com