

# **ELECTRIC OBSTERIC BED HBD35-EO**



# ELECTRIC OBSTERIC BED HBD35-EO

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## HBD35-EO ELECTRIC OBSTERIC BED

Imported LINAK motor various actions are controlled by remote controller

Super mute big foot wheel which is easily moved and equipped with center control breaking system.

Equipped with stainless steel sewage basin

The mattress adopts high density sponge and manufactured by high quality artificial leather without seam

With various complete accessories

Standard Accessories: A pair of leg support frame, a stainless steel sewage basin, a stand for transfusion, a remote control, a power wire



## SPECIFICATIONS

Model	HBD35-EO
Length	1950+50 mm
Width	800+50 mm
Height Adjustment	(550~800)+50 mm
Trendelenburg/rev-trend	$\geq 5^\circ / \geq 12^\circ$
Backboard upward	$\geq 60^\circ$
Foot board upward	$\geq 40^\circ$
Foot board downward	$\geq 90^\circ$
Standard accessory	A pair of leg support, a sewage basin, a power cord, a IV pole, a remote controller
Input power	150+10%
Power Supply	220+22V, 50+1Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)