

MANUAL MULTI- FUNCTION GYNEACOLOGICAL TABLE HBD36-MM



MANUAL MULTI-FUNCTION GYNEACOLOGICAL TABLE HBD36-MM

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gyneacological abdominal operation, Vulva operation.

Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

HBD36-MM MANUAL MULTI-FUNCTION GYNEACOLOGICAL TABLE

Material: Steel Powder Coated S.S, Color: Blue



SPECIFICATIONS

| Model | HBD36-MM |
|--------------------------|-----------------------|
| Function | Multi |
| Back-rest lifting angle | $\geq 20^\circ$ |
| Back section up adjusted | $\geq 60^\circ$ |
| Seat board upwards | $\geq 8^\circ$ |
| Seat board downwards | $\geq 10^\circ$ |
| Load Bearing | 250 kg |
| Height adjustment | ≤ 650 mm |
| Tabletop Dimension | 1240Lx600W mm |
| Upright Dimension | 500Lx600W mm |
| Dimension | 1750Lx600Wx650~900 mm |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com