

# ELECTRIC GYNECOLOGICAL TABLE HBD38-1270



# ELECTRIC GYNECOLOGICAL TABLE HBD38-1270

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## HBD38-1270 ELECTRIC GYNECOLOGICAL TABLE

Various actions are controlled by Imported LINAK motor

Hand-held remote control and pedal remote control, independent of each other, easy to operate it

With concealed stainless steel sewage waste basin

The mattress is wrapped with imported artificial leather. High density memory sponge is formed in one piece, with no seam

Equipped with brake device which is easy for moving the table

The base is made up of high strength ABS, which is corrosion resistant and easy to clean

Leg plate can be lifted electrically

Concealed stainless steel auxiliary table

With function of auto restoration



## SPECIFICATIONS

Model	HBD38-1270
Table Length	1270+50 mm
Table width	550+30 mm
Height Adjustment	(700-1000)+50 mm
Trendelenburg/rev-trend	$\geq 5^\circ / \geq 12^\circ$
Backboard downward	$\geq 40^\circ$
leg board reflex	$\geq 30^\circ$
Working area	(350x480)+20 mm
Standard accessory	A pair of arc shaped handle, a pair of leg support, a stainless steel sewage basin, a power cord, a hand held remote controller, a foot controller
Main Voltage AC	220V+22V, 50Hz+1Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)