

# **ELECTRIC GYNECOLOGICAL TABLE HBD38-1300**



# ELECTRIC GYNECOLOGICAL TABLE HBD38-1300

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## HBD38-1300 ELECTRIC GYNECOLOGICAL TABLE

Adopt high density memory sponge and is once formed, with no seam  
Hydraulic lifting by foot pedal, other actions are controlled by pneumatic  
Imported Y shaped sealing ring and durable  
Stainless steel concealed auxiliary table  
The base is made up of high strength ABS  
Equipped with foot brake device which is easy for moving the table



### SPECIFICATIONS

Model	HBD38-1300
Table Length	1300+50 mm
Table width	550+30 mm
Height Adjustment	(700 950)+50 mm
Trendelenburg/rev-trend	$\geq 3^\circ / \geq 10^\circ$
Dimension of auxiliary platform	(350x450)+20 mm
Backboard downward	$\geq 70^\circ$
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)