

# ELECTRIC GYNECOLOGICAL TABLE HBD38-1450



# ELECTRIC GYNECOLOGICAL TABLE HBD38-1450

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## HBD38-1450 ELECTRIC GYNECOLOGICAL TABLE

Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with hidden sewage basin

With foot control, can adjust the movements easily

Imported LINAK gear motor, without noise, have a long service life



### SPECIFICATIONS

Model	HBD38-1450
Table Length	1450+50 mm
Table width	520+30 mm
Height Adjustment	(400 900)+50 mm
Trendelenburg/rev-trend	$\geq 15^{\circ}/\geq 25^{\circ}$
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Power Supply	220V 50Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)