

ELECTRIC GYNECOLOGICAL TABLE HBD38-2000



ELECTRIC GYNECOLOGICAL TABLE HBD38-2000

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

HBD38-2000 ELECTRIC GYNECOLOGICAL TABLE

Imported LINAK motor various actions are controlled by foot controller

The base seat and upright column covering are made of high quality 304 stainless steel which is anti corrosive and can be easily clean

External folding auxiliary table is for easy dismantling

Equipped with stainless steel sewage basin and other accessories

Foot treadle brake device which is easy for moving obstetric table



SPECIFICATIONS

Model	HBD38-2000
Table Length	2000+50 mm
Table width	600+10 mm
Height Adjustment	(650-900)+50 mm
Trendelenburg/rev-trend	$\geq 10^{\circ}/\geq 12^{\circ}$
Backboard upward	$\geq 50^{\circ}$
Leg board outward	$\geq 90^{\circ}$
Working area length/width	600x600 mm
Standard accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of foot arm rods, a stainless steel sewage basin, a foot controller

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com