

HYDRAULIC TRANSFER STRETCHER HBD81-HS



HYDRAULIC TRANSFER STRETCHER HBD8I-HS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8I-HS HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher

Five functions available by pressing the foot-pedals

Right arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easily control the stretcher direction

Standard accessories: Mattress, I.V.Pole



SPECIFICATIONS

Model	HBD8I-HS
Back-rest	0 – 85°
Height range (bed board to ground)	813 – 590 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Trendelenburg	0 – 12°
Reverse trendelenburg	0 – 12°
Overall dimension	1906x743 mm
Net Weight	82 kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com