

# HYDRAULIC TRANSFER STRETCHER HBD8J- HS1



# HYDRAULIC TRANSFER STRETCHER HBD8J-HS1

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## HBD8J-HS1 HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher, 5 functions available by hydraulic and gas spring and crank system together

At the head-end, there is a pair of push- handles, can be folded up and down

Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily, Bed board for X-Ray Examination



## SPECIFICATIONS

Model	HBD8J-HS1
Back-rest	0 – 85°
Height range (bed board to ground)	860 – 630 mm
Knee Gatch	0°- 50°
Trendelenburg	0 – 13°
Reverse trendelenburg	0 – 13°
Standard accessory	IV pole
Bed Board	1932x630 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	2050x745 mm
Net Weight	105 kg
Height adjustment	835-600 mm

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)