

MANUAL WHEELCHAIR HBD91- MW5



MANUAL WHEELCHAIR HBD91-MW5

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

HBD91-MW5 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



SPECIFICATIONS

Model	HBD91-MW5
Wheelchair length	1250 mm
Opened wheelchair width	690 mm
Wheelchair Height	1210 mm
Closed wheelchair width	365 mm
Seat Width	450 mm
Front wheel	8" inch
Rear wheel	24" inch
Seat Height	520 mm
Seat depth	455 mm
Backrest Height	450 mm
Armrest Height	240 mm
Distance between armrest	450 mm
Height between footrests and ground	100 mm
Net Weight	25.4 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com