

# ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR HBD93-123A



# ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

## HBD93-123A

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient transfer, patient lift, toilet-bath chair, Patient Indoor-outdoor locomotion.

Also known as Motorized wheelchair, Electric-powered wheelchair.

## HBD93-123A ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

This is a kind of patient lift

it is electric lift

height of sling can be adjustable



### SPECIFICATIONS

Model	HBD93-123A
Material	Iron, Plastic, Nylon
Maximum loading weight	150 kg
Lifting Range	25 cm, from 40 cm to 65 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 55 cm
It suites for bed	Height of bed from 11 cm to 63 cm
Waterproof level	IP44
Dimensions	123x72.5x54.5 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, rechargeable

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)