

PATHOLOGICAL SLIDE SCANNER HTP51-100



PATHOLOGICAL SLIDE SCANNER HTP51-100


Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision.

Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens.

Also known as Slide printer, Cassettes printer.

HTP51-100 PATHOLOGICAL SLIDE SCANNER

The PrimeHisto XE is a histological slide scanner. With easy manipulation/setting, PrimeHisto XE allows the histological samples to be digitized in a matter of seconds. Ability to be suitable for most histological slide samples and can scan the entire histological slide image in a scanning area of 24.3x36.5mm. Also, has a high resolution of up to 10,000 DPI.



Product Image Coming Soon

SPECIFICATIONS

Model	HTP51-100
Resolution	10000x10000 DPI
Dynamic range	3.8 OD
Scan area	24.3x36.5 mm
Scanning medium	75x25 mm section
Photo formats	TIFF, JPEG
Light	transmits LED white light
Operating systems	Windows XP/Vista / 7/8/10, Mac OS X 10.5+
Volume	260x155x70 mm
Weight	1.05 kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com