

MINI COVER SLIPER HTP51-15



MINI COVER SLIPER HTP51-15

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision. Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens. Also known as Slide printer, Cassettes printer.

HTP51-15 MINI COVER SLIPER

Covering machine is a fast curing device for sealed glass slides. This machine adopts advanced microcomputer technology, which makes the equipment more accurate, stable and anti-interference ability.



Product Image Coming Soon

SPECIFICATIONS

Model	HTP51-15
Slide seat	15
Speed	30 tablets/min
Rated voltage	AV 220 V \pm 10 %
Power frequency	50 \pm 1 Hz
Rated power	\leq 100W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@m.labstac.com | Website: labstac.com