

# SLIDE PRINTER HTP51-150




# SLIDE PRINTER HTP51-150

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision. Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens. Also known as Slide printer, Cassettes printer.

## HTP51-150 SLIDE PRINTER

Laser printing and easy for maintaining. Capability of linking with the management system of the pathology department. Print online or off-line. The writing is clear and delicate, never wears out or fades, which meets the characteristic requirements of keeping slides for long-term reference. No need for preheating, can be immediately used after starting up, safe and noiseless, and free of Maintenance.



### SPECIFICATIONS

Model	HTP51-150
Capacity	150 Slides
Average printing speed	≤ 3S / slide
Printing accuracy	≥ 1250 pi
Maximum Capacity	≥ 150 slides
Minimum Capacity	≥ 7 slides
Ambient Humidity	45 – 85 %
Dimensions	465x200x361 mm (LengthxWidthxHeight)
Laser Power	Pulse laser ≤ 2 W
System interface	Integrate with LIS/HIS/P IS
Printing Methods	Print online or off-line
Loading method	Cartridge clip type
External device	Support Bar-code scanner inputting
Consumable items	5×75 mm, 1 - 1.2 mm glass slide

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)