

# CASSETTE PRINTER HTP51-75



# CASSETTE PRINTER HTP51-75

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision.

Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens.

Also known as Slide printer, Cassettes printer.

## HTP51-75 CASSETTE PRINTER

Laser printing and easy to maintain. Capability of linking with the management system of the pathology department. Laser etching on the surface of the embedding cassettes, no need for consumables. No need for preheating, can be immediately used after starting up, safe and noiseless, free of

Maintenance.



### SPECIFICATIONS

Model	HTP51-75
Printing method	Cassette Laser Printing
Average printing speed	$\leq 3S$ / Cassette
Printing accuracy	$\geq 1250$ pi
Maximum Capacity	$\geq 75$ Cassettes
Minimum Capacity	$\geq 1$ slides
Power	100-240 ACV 50/60 Hz
Standby power	50 W
Rated Power	100 W
Dimensions	395x280x330 mm(LengthxWidthxHeight)
Laser Power	Pulse laser $\leq 2$ W
System interface	Integrate with LIS/HIS/P IS
Printing Methods	Print online or off-line
Adhesion testing	Pull 1x10x30 mm 3M tape 3 times without fade
External device	Support Bar-code scanner inputting
Consumable items	Standard cassette

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)