

BULLET ICE MAKER ICM11-020



BULLET ICE MAKER ICM11-020

Compact machine with ice making capacity and innovative technology. It offers durability, efficiency and saving. Using quality assured components, it gives excellent functionality and low maintenance.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

ICM11-020 BULLET ICE MAKER

Vertical stirring engine model

Refrigeration system produces uniform shape ice bullets

Bell-shaped ice with a hexagon hole

Highly durable stainless steel 304 2B built-up

Automatic controller increases working reliability

Large ice storage bin

Eco-friendly, energy saving non-Freon refrigerant



SPECIFICATIONS

Model	ICM11-020
Capacity	20 kg/24 h
Storage Bin Capacity	6 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/32 kg
Power	220 W
Power Supply	220V /110 V, 50 / 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2600505005	Water filter

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com