

ICE FLAKER ICM31-020



ICE FLAKER ICM31-020

Suits a variety of applications including fisheries, laboratories and industry. Once formed, Ice flakes are dry, clean, beautiful in shape. Relatively big contact area and good mobility can ensure their full contact with refrigerated objects to realize good preservation effect. Great adaptability, sanitary and convenient.

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

ICM31-020 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



SPECIFICATIONS

Model	ICM31-020
Capacity	20 kg/24 h
Storage Bin Capacity	10 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/31 kg
Power	280 W
Power Supply	220 V, 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2600505005	Water filter

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com