

GENERAL PURPOSE INCUBATOR INC073-016



GENERAL PURPOSE INCUBATOR INC073-016

INC073-016 GENERAL PURPOSE INCUBATOR

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control

Air circulation system with specific air flow channel ensures a good temperature uniformity performance

Silicon door gasket with long lifetime, and easy to change

A damper adjustment in the front ensures the gas convection enough in working chamber



SPECIFICATIONS

Model	INC073-016
Capacity	16 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	850 W
Internal Dimension	250Wx260Dx250H mm
External Dimension	530Wx370Dx420H mm
Shelves	2 pcs
Timing Range	0~9999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com