

GENERAL PURPOSE INCUBATOR INC073-136



GENERAL PURPOSE INCUBATOR INC073-136

INC073-136 GENERAL PURPOSE INCUBATOR

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control

Air circulation system with specific air flow channel ensures a good temperature uniformity performance

Silicon door gasket with long lifetime, and easy to change

A damper adjustment in the front ensures the gas convection enough in working chamber



SPECIFICATIONS

Model	INC073-136
Capacity	136 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	2050 W
Internal Dimension	550Wx450Dx550H mm
External Dimension	840Wx580Dx730H mm
Shelves	2 pcs
Timing Range	0~9999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com