

# GENERAL PURPOSE INCUBATOR INC073-220



# GENERAL PURPOSE INCUBATOR INC073-220

## INC073-220 GENERAL PURPOSE INCUBATOR

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control

Air circulation system with specific air flow channel ensures a good temperature uniformity performance

Silicon door gasket with long lifetime, and easy to change

A damper adjustment in the front ensures the gas convection enough in working chamber



## SPECIFICATIONS

Model	INC073-220
Capacity	220 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	2450 W
Internal Dimension	600Wx500Dx750H mm
External Dimension	880Wx630Dx930H mm
Shelves	2 pcs
Timing Range	0~9999 min

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)