

GENERAL PURPOSE INCUBATOR INC073-620



GENERAL PURPOSE INCUBATOR INC073-620

INC073-620 GENERAL PURPOSE INCUBATOR

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control

Air circulation system with specific air flow channel ensures a good temperature uniformity performance

Silicon door gasket with long lifetime, and easy to change

A damper adjustment in the front ensures the gas convection enough in working chamber



SPECIFICATIONS

Model	INC073-620
Capacity	620 L
Electrical requirement	380 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	4000 W
Internal Dimension	840Wx600Dx1355H mm
External Dimension	980Wx800Dx1880H mm
Shelves	4 pcs
Timing Range	0~9999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com