

GENERAL PURPOSE INCUBATOR INC080-016



GENERAL PURPOSE INCUBATOR INC080-016

INC080-016 GENERAL PURPOSE INCUBATOR

Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	INC080-016
Capacity	16 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	200 W
Interior Dimension	250Wx260Dx250H mm
Exterior Dimension	530Wx480Dx420H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com