

GENERAL PURPOSE INCUBATOR INC080-035



GENERAL PURPOSE INCUBATOR INC080-035

INC080-035 GENERAL PURPOSE INCUBATOR

Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	INC080-035
Capacity	35 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	200 W
Interior Dimension	340Wx320Dx320H mm
Exterior Dimension	620Wx490Dx490H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com