

# GENERAL PURPOSE INCUBATOR INC080-160



# GENERAL PURPOSE INCUBATOR INC080-160

## INC080-160 GENERAL PURPOSE INCUBATOR

Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature



## SPECIFICATIONS

Model	INC080-160
Capacity	160 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	600 W
Interior Dimension	500Wx500Dx650H mm
Exterior Dimension	790Wx630Dx810H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)